



Charisma® Topaz

Innovative chemistry for day-to-day natural restorations.

Giving a hand to oral health.





Savour the effortlessness of easy handling.

- The soft, yet low-sticky consistency of Charisma Topaz improves the efficiency during modelling thanks to its increased stability.
- The extended working time even under operating and ambient light allows for easy, stress-free modelling in any indication.
- Our filler system with nano particles facilitates polishing and ensures both a high, long-lasting gloss and high patient satisfaction.



Benefit from superior physical properties.

- We can allow ourselves to promise reliability and quality with Charisma Topaz, which are the keyfactors for long-lasting restorations. Its patented TCD matrix combines high strength with minimal shrinkage, both are important material properties which can minimise the risk of fractures and secondary caries.
- Its good radiopacity enables a safe and reliable diagnosis.
- You do not need to change any of your proceeding habits, as Charisma Topaz is compatible to all conventional bonding materials and techniques.



Take advantage from versatility in layering.

- Use the easy single-shade layering with universal shades for your less complex cases and reduce your layering efforts to a minimum.
- The multi-shade layering possibilities allow natural layering also of complex cases. Choose from three levels of opacity plus bleach shades.
- With Charisma Topaz, restorative margins are virtually undetectable due to the excellent shade match.



The patented TCD-matrix is free of BPA related monomers Further information kulzer.com/charisma-material








Charisma® Topaz

Innovative chemistry for day-to-day natural restorations.


Easy handling meets innovative formula for aesthetic versatility.



Charisma® Topaz

Down-to-earth natural beauty.

Savour the effortlessness of easy handling.



Five minutes after a diamond form was pressed into the composite, Charisma Topaz still shows sharp and stable impressions whereas other composites show unfavorable slumping.

Source: Kulzer CRE, Hanau

Benefit from superior physical properties.



Source: Kulzer R&D, Wehrheim



Pictures by courtesy of Prof. Dr. med. dent. Claus-Peter Ernst, Mainz, Germany

Charisma® Topaz

The three-in-one solution.

Take advantage from versatility in layering.



Before



Replacement of an gold inlay using a single-shade layering technique



Before



Multi-shade layering technique in the aesthetic zone

Charisma Topaz product range

Charisma Topaz Syringe Basic Kit + GLUMA Bond Universal

1 x 4 g syringe each of the shades A2; A3; A3.5; B2; CL; OM 1 x 4 ml bottle of GLUMA Bond Universal; 1 x 2.5 ml syringe of GLUMA Etch 35 Gel; Layering Guide, Pictorial Cards Art.-No.: 66067179

Art.-No.: 66067180

Charisma Topaz Syringe Master Kit + GLUMA Bond Universal

1 x 4 g syringe each of the shades A1; A2; A3; A3.5; B2; C2; CL; OL; OM; OD 1 x 4 ml bottle of GLUMA Bond Universal; 1 x 2.5 ml syringe of GLUMA Etch 35 Gel; Shade Guide, Layering Guide, Pictorial Cards



Charisma Topaz Shade Guide Art.-No.: 66067181





Shade	Syringe
	1 x 4 g
Universal shades	
A1	66067163
A2	66067164
А3	66067165
A3.5	66067166
A4	66067167
B1	66067168
B2	66067169
В3	66067170
C2	66067171
Opaque shades	
OL (Opaque Light)	66067172
OM (Opaque Medium)	66067173
OD (Opaque Dark)	66067174
Incisal shades	
CL (Clear)	66067175
CO (Clear Opal)	66067176
Bleach shades	
BL (Bleach Light)	66067177
OB (Opaque Bleach)	66067178

The scope of bonding in one bottle, that is GLUMA® Bond Universal.

You do not want to change your bonding habits? You would like to keep your bonding routines as easy as possible? Based on our long-term bonding expertise, we developed GLUMA Bond Universal – a universal bonding answering to all your bonding needs. It is a reliable and highly effective adhesive that ensures successful long-term restorations.

One bottle of our genuinely universal adhesive brings you manifold benefits:

- One bottle for any bonding technique: total-etch, selective etch, and self-etch.
- One bottle for all dental materials.
- One bottle for a broad range of indications of direct and indirect restorations and intraoral repairs.
- One bottle for easy and safe handling with state-of-the-art performance.





Contact in Germany

Kulzer GmbH Leipziger Straße 2 63450 Hanau, Germany info.dent@kulzer-dental.com

Filtek® Z250, Filtek® Z350 and Filtek® Supreme XTE are trademarks of 3M Company • Esthetx® HD is a trademark of Dentsply Sirona Inc. • Herculite® XRV Ultra™ is a trademark of Kerr Corporation • Gradia® Direct, G-aenial® Posterior and G-aenial® Anterior are trademarks of GC Corporation • Tetric® N-Ceram is a trademark of Ivoclar Vivadent AG • Estellte® Sigma Quick is a trademark of Tokuyama Dental Corporation Inc.