



GLUMA® Self Etch

Strong bonding for life.

One step ahead.

Based on our long-term bonding expertise, we have developed the self-adhesive system GLUMA Self Etch. This all-in-one bonding agent etches, primes, bonds and desensitises in one single step. Thus, it meets your everyday requirements and demonstrates all the characteristics which are important for an uncomplicated and state-of-the-art bonding agent.

With this self-adhesive one-bottle agent you benefit from:

- » The fast and effective one-coat application for etching, priming, bonding plus desensitising.
- » Excellent bond strength and marginal sealing.
- » Optimal usage in hypersensitive areas.
- » Reliable, scientifically proven performance.



Giving a hand to oral health.



KULZER
MITSUI CHEMICALS GROUP

Filename: KULZ1708_11975_GLUMA_Self_Etch_Product_Sheet_EN_Kulzer	Kulzer Toolbox No: W19087	Size closed (EU, USA): EU 210x297mm	Number of pages 2 pages
SAP No.: -	ORT-Version: 02	Speciality in printing/Spotcolours: -	Lacquer: -
Marketing Manager: S. Oyegbesan	Producer: ORT	printing profile: iso coated v2 (eci)	Date 09-12-2017
			Kulzer final release 09-12-2017

GLUMA® Self Etch



Etching, priming and bonding in only one step with a desensitising effect!

GLUMA Self Etch is applied in only one step. A shiny surface is easily achieved by air drying and provides the base for outstanding, durable restorations. The excellent stability of GLUMA Self Etch allows storage at room temperature.

Indications for GLUMA Self Etch:

- Bonding of direct light-cured composite restorations (including compomers).
- Bonding indirect restorations in combination with a light-curing luting cement: porcelain and composite restorations (inlays, onlays, veneers, crowns).
- Treatment of hypersensitive tooth regions.

Application Step by Step:



Source: Degrange M, University of Paris, France 2007. Documentation available. GLUMA Self Etch is equivalent to the tested product.

Art. code	Article
66046246	GLUMA Self Etch Bottle Refill 1 x 4 ml Bottle, 50 Application Tips, Mixing Well, Pictorial Card
66046247	GLUMA Self Etch Bottle Value Pack 3 x 4 ml Bottle Pictorial Card



OptiBond® is a trademark of Kerr Corporation • XP Bond® is a trademark of Dentsply De Trey GmbH • Scotchbond® 1XT is a trademark of 3M Company • Admira® is a trademark of Voco GmbH • Clearfil® SE Bond and Clearfil S3 Bond® are trademarks of Kuraray Co. Ltd. • AdheSE® is a trademark of Ivoclar Vivadent, Inc. • G-Bond® is a trademark of GC America

© 2017 Kulzer GmbH. All Rights Reserved.

Kulzer GmbH · Leipziger Straße 2 · 63450 Hanau · Germany · kulzer.com

W19087 EN 07/2017 ORT